



NEW
Fits for all nitrous
oxide systems

NITROUS OXIDE / OXYGEN SCAVENGING SYSTEM



Double Nasal Hoods and Scavenging Circuits

Bringing Nitrous Oxide / Oxygen Sedation
to the next level with multiple innovations

Baldus
Sedation

Nitrous Oxide / Oxygen Sedation



360° SERVICE

- » Premium Manufacturer
- » Customer Service
- » Marketing Power Support
- » The only German manufacturer of nitrous oxide sedation systems for dentists
- » Market leader in Germany
- » Global partners in over 50 countries
- » Used by medical professionals in hospitals, surgical clinics, laboratories and dental practices

Different at First Sight



Adult



Pediatric



PROFILE

Lower profile for a better access to the patient



MATERIAL


Softer and more flexible for an exceptional seal and comfort



DESIGN

Oval design of the exhalation valve for ideal scavenging, US patented*

Innovations at a Glance



MATERIAL



- » Soft inner hood improves seal and patient comfort
- » Clear outer hood sits closer to the patient's face, increasing scavenging efficiency of exhaled gases and allowing visual monitoring of patient's breathing
- » Latex-free material protects patients and staff who have latex allergies

COMPATIBLE

Universal fit with other double scavenging circuits allows for use of your existing systems with the new Baldus Sedation Double Nasal Hoods




DESIGN

- » Double hood design minimizes ambient nitrous oxide in operatory
- » Baldus Sedation Double Nasal Hood utilize an oval exhalation valve that directs exhaled gases to the inner hood's built-in scavenger "flow channels" maximizing efficiency


PROFILE



- » **Accurate fit!** The inner hood seals to the face to ensure proper gas delivery to the patient
- » **Size Matters!** Smaller profile provides improved "sight line" visibility for the dental staff and therefore allows better access to the oral cavity



Baldus Sedation Double Nasal Hood compared to other available double twice nasal hoods



.....Peer 1
.....Peer 2

BALDUS Sedation Double Nasal Hoods

- » Colorful and lightly scented single use hoods engage and delight patients of all ages
- » Pediatric and adult sizes available in a variety of scents

12 HOODS INDIVIDUALLY WRAPPED IN A BOX

Part no.	Scent	Product description, contents and price
Adult		
16030-ABlue	Blueberry Dance	Double Nasal Hood Single Use, 12 Hoods individually wrapped in a box € 99,00
16030-AE	Strawberry Delight	
16030-ABub	Bouncy Bubble	
16030-AA	Happy Apple	
16030-AV	Chilla Vanilla	
16030-AG	Pure Ghost	
Pediatric		
16030-PBlue	Blueberry Dance	Double Nasal Hood Single Use, 12 Hoods individually wrapped in a box € 99,00
16030-PE	Strawberry Delight	
16030-PBub	Bouncy Bubble	
16030-PA	Happy Apple	
16030-PV	Chilla Vanilla	
16030-PG	Pure Ghost	



Autoclavable hoods also available

Nitrous Oxide / Oxygen Scavenging Circuits



- » Fully autoclavable scavenging circuit provides for safer delivery and monitoring

DOUBLE HOOD SCAVENGING CIRCUIT WITH SINGLE USE NASAL HOODS STARTER KIT

Part no.	Product description, contents and price
16055	Double Hood Scavenging Circuit including Single Use Nasal Hoods, 3 Adult (Scent Blueberry Dance) and 3 Pediatric (Scent Strawberry Delight) € 449,00
16056	Double Hood Scavenging Circuit w/ Vacuum Control Block including Single Use Nasal Hoods, 3 Adult (Scent Blueberry Dance) and 3 Pediatric (Scent Strawberry Delight) € 599,00



MANUFACTURER
Baldus Sedation GmbH & Co. KG
In der Langfuhr 32
56170 Bendorf, Germany

Phone: +49 261 9638926-66
Fax: +49 261 9638926-20
E-mail: info@baldus-sedation.de
Website: www.baldus-sedation.de/en/

Baldus
Sedation